



MATERIAL MEDICO Y QUIRURGICO
EXCLUSIVAS MÉDICAS
informacion@exclusivasmedicas.com
Tlf.: 00/34/ 986 508 028 Fax: 00/34/ 986 508 027

www.carewell.com.cn



CUS-9618D

Digital Ultrasound Scanner

More Than Excellent Image



CUS-9618D

Digital Ultrasound Scanner



► Main features:

- Full digital beam-forming and dynamic aperture technology ensure excellent image quality
- 96 element, 32 channels, dual sockets
- scanning angle 8-16 level adjustable
- gamma correction, 8-segment GTC
- automatic report generation(normal/OB)
- abundant software packages for full clinic applications such as abdominal, obstetric, gynecology, cardiology, urology, small parts, etc.
- ergonomic design, user-friendly

► Specifications:

Scanning: electronic convex, linear and micro-convex
Display: B, B/B, B/M, 4B
Image conversion: up/down, right/left, black/white
Notation: real-time clock, date, ID, gender, age, focus, hospital, doctor
Scanning depth: 220mm(max)
Gray scale: 256
Display: 10" SVGA
Video output: PAL-D
Body marks: 15
Cine-loop: 829 frames(max)
Permanent image storage: 32 frames
Image magnification: x1.0, x1.2, x1.5, x2.0
Picture in picture: 1:2
Pseudo-color: 0-7
Measurement of distance, circumference, area, volume, HR, EF, EDD
USB2.0 interface

► Standard probe

C3-1/R60/3.5MHz multi-frequency convex



R60/3.5MHz Convex

► Other options

Video printer, laser printer, trolley cart, acupuncture bracket & needle



► Option probe

L3-1/7.5MHz multi-frequency linear
EL3-1/7.5MHz multi-frequency endorectal linear
EC1-1/R13/6.5MHz multi-frequency transvaginal convex
C1-6/20R/5.0MHz micro-convex



R20/5.0MHz Micro-convex L40/7.5MHz Linear R13/6.5 Transvaginal L40/7.5 Endorectal

Specifications subject to change without prior notice



MATERIAL MEDICO Y QUIRURGICO
EXCLUSIVAS MÉDICAS
informacion@exclusivasmedicas.com
Tlf.: 00/34/ 986 508 028 Fax: 00/34/ 986 508 027

